

NEOTEMP

**Precision Neonatal Temperature Monitoring
for GE/Ohmeda Medical and Dräger Medical/Air Shields
Incubators & Infant Radiant Warmers.**









- Temp Accuracy** Baby-sensitive thermistor sensor guarantees temperature reading sensitivity with a stable, highly accurate measurement.
- Compatibility** Choose from a comprehensive range of sensors designed for GE / Ohmeda Medical and Dräger Medical / Air Shields incubators and radiant warmers.
- Ease of Use** Molded plug assembly provides direct connection to the incubator/ warmer eliminating the need for an interface cable.
- Sensor Safety** Flexible lead wire facilitates sensor placement directly on the infant further ensuring precise temperature reading.
- Gentle** Ensures patient safety with a one inch round baby sensitive hypoallergenic adhesive cover. A reflective Mylar adhesive cover provides insulation from heat and external sources.

NEONATAL Replacement Adhesive Mylar Disc Covers		
10 - 1900 - 300	Mylar Adhesive Replacement Disc Covers -	300 / Box

NEOTEMP

Disposable Neonatal Skin Temperature Sensors

(Packaged 20/box, 60/case)

		Part #	Description	Substitute For
GE/Ohmeda Medical/Ohio		10 - 1930 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for GE Care Plus, I.C., CG Incubators & Infant Warming Systems (IWS) 3000, 3300, 3500, 5000 & Panda Warmers	Replaces OEM # LA003, LA005, # 6600 - 0196-700 & # 6600 - 0208-700
		10 - 1930 - 060 Packaged 60/case		
GE/Ohmeda Medical/Ohio		10 - 1950 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for GE GIRAFFE OmniBed, Giraffe Incubators & Warmers, New Panda Warmers	Replaces OEM # 6600 - 0873-700 & # 6600 - 0874-700
		10 - 1950 - 060 Packaged 60/case		
Dräger Medical/Air-Shields		10 - 1910 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for ISOLETTE Incubators (Models C-100, C-450, C-550) & PM-78 Radiant Warmers	Replaces OEM #68-209-20/MU06933 & #68-209-30/MU06935
		10 - 1910 - 060 Packaged 60/case		
		10 - 1960 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for ISOLETTE C2000 & 8000 Incubators, VERSALET 7700 Care Center	Replaces OEM # 83-104-15/MU12551
		10 - 1960 - 060 Packaged 60/case		
	10 - 1965 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for CALEO Incubators & Babytherm Warmer 8004/8010	Replaces OEM # MX11000	
	10 - 1965 - 060 Packaged 60/case			
	10 - 1970 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for RESUSCITAIRE Radiant Warmers	Replaces OEM # 81-300-006, # 81-300-008/MU11228 & # 81-300-009/MU11229	
	10 - 1970 - 060 Packaged 60/case			
Atom Medical		10 - 2000 - 020 Packaged 20/box	NOVATEMP Skin Temperature Sensor For Atom Medical Incubators & Infant Warmers	Replaces OEM #60882, #60883, & #60084 w / connector