

NeoTemp Neonatal Skin Temperature Sensors

Choose from a complete range of NEOTEMP sensors compatible with Dräger Medical/Air Shields & GE /Ohmeda Medical incubators and infant radiant warmers.

		Part #	Description	Substitute For
GE/Ohmeda Medical/Ohio		10 - 1930 - 020 Packaged 20/box 10 - 1930 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for GE Care Plus, I.C., CG Incubators & Infant Warming Systems (IWS) 3000, 3300, 3500, 5000 & Panda Warmers	Replaces OEM # LA003, LA005, # 6600 - 0196-700 & # 6600 - 0208-700
		10 - 1950 - 020 Packaged 20/box 10 - 1950 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for GE Omni-Bed GIRAFFE, GIRAFFE Incubators & New PANDA Warmers	Replaces OEM # 6600 - 0873-700 & # 6600 - 0874-700
Dräger Medical/Air-Shields		10 - 1910 - 020 Packaged 20/box 10 - 1910 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for DM ISOLETTE Incubators (Models C-100, C-450, C-550) & PM-78 Radiant Warmers	Replaces OEM #68-209-20/MU06933 & #68-209-30/#MU06935
		10 - 1960 - 020 Packaged 20/box 10 - 1960 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for DM ISOLETTE C2000 Incubators & Versalet 7700 Care Center	Replaces OEM # 83-104-15/#MU12551
		10 - 1965 - 020 Packaged 20/box 10 - 1965 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for DM Caleo Incubators	Replaces OEM # MX11000
		10 - 1970 - 020 Packaged 20/box 10 - 1970 - 060 Packaged 60/case	NEOTEMP Neonatal Skin Temperature Sensor for Resuscitaire Infant Radiant Warmers	Replaces OEM # 81-300 - 006, # 81-300 - 008/MU11228 & # 81-300 - 009 / MU11229
		10 - 1600 - 030 Packaged 30/box	NovaTemp Skin Temperature Sensor for AirBorne Transport and Dräger Transport Incubators	400 YSI Connector Requires 400 Series Cable #11-CA-402-12

NOVAMED USA • The Specialist in Neonatal Temperature Monitoring
800 - 425-3535/914-789-2100 • www.NovaMed-USA.com