



PREMIER Fiber Optic Laryngoscopes

Integrated Fiber Optics for Enhanced
Light Transmission, Improved Illumination



Reusable Integrated Fiber Optics
Xenon Technology – Green System



Advanced Engineering for Maximum Light Transmission NOVALITE PREMIER Fiber Optic Laryngoscopes

Reusable Integrated Fiber Optics for Advanced Performance

Designed with high quality stainless steel construction, enlarged fiber optic bundle for improved illumination and safer intubations.

NOVALITE PREMIER Fiber Optic Laryngoscopes utilize the universal Green System to ensure interchangeability with competitive brands.

- PREMIER MacIntosh blades feature a horizontal 5mm diameter fiber optic bundle with 7000+ individual micro-fibers
- XENON technology provides focused light transmission [+40% compared to halogen bulbs].
- Smooth, integrated design facilitates ease of sterilization, helps avoid cross-infection.
- Low profile design ensures safer intubations, enhanced patient comfort.
- High quality Stainless Steel construction ensures long life and durability.
- Available in a complete selection of MacIntosh and Miller blades.
- Compatible with all ISO 7376 Green System handles.



PREMIER blades feature low profile smooth lines for safer intubations.

Integrated Fiber Optics



PREMIER blade design features smooth lines from end to end, eliminating cracks and crevices for ease of maintenance and longer life.

Safety



Integrated fiber optic bundle promotes ease of cleaning and sterilization, eliminating risk of cross-infection. This is important to enhance patient safety and improve risk management.

Green System



Utilizes the universal green system to ensure compatibility with all green system handles.



PREMIER Fiber Optic Integrated Laryngoscope Blades

Reusable Stainless Steel Green System, Satin Finish



MACINTOSH PREMIER Integrated Fiber Optic Reusable Blades	
Size / Description	Order #
0, Neonate	17 - PIB - MAC - 0
1, Infant	17 - PIB - MAC - 1
2, Child	17 - PIB - MAC - 2
3, Medium Adult	17 - PIB - MAC - 3
4, Adult	17 - PIB - MAC - 4
5, Large Adult	17 - PIB - MAC - 5



MILLER PREMIER Integrated Fiber Optic Reusable Blades	
Size / Description	Order #
00, Premie	17 - PIB - MIL - 00
0, Neonate	17 - PIB - MIL - 0
1, Infant	17 - PIB - MIL - 1
2, Child	17 - PIB - MIL - 2
3, Medium Adult	17 - PIB - MIL - 3
4, Adult	17 - PIB - MIL - 4





Flip Tip Fiber Optic MacIntosh Blades, Reusable

Ideal for difficult intubations with flexible tip adjustable through 70°

Size / Description	Order #
2, Child	14 - 200FOFT - MAC - 2
3, Medium Adult	14 - 200FOFT - MAC - 3
4, Adult	14 - 200FOFT - MAC - 4



Fiber Optic Green System Reusable Laryngoscope Handles		
Size	Batteries	Order #
Medium/Adult	2C	14 - 600FO - 01
Penlight/Pediatric	2AA	14 - 600FO - 02
Stubby	2AA	14 - 600FO - 03



Xenon Light Technology Replacement Bulbs	
Description	Order #
Ultra Brite Xenon Bulbs 2.7V (6 pack)	14 - 700FO - 02 - 06